

Development of a novel retinal tamponade to replace gas and oil

Tomasz P Stryjewski MD
Tallman Eye Associates
Boston, Massachusetts

Retina Society 2020

Financial Disclosure



Co-founder and equity co-owner with Tony Stefater MD PhD in a start-up company (Pykus Therapeutics, Inc., Cambridge, MA) which is developing the presented technology for clinical use



Summary



- Gas/oil tamponade is a considerable burden for patients
- Many promising alternatives have not advanced to the clinic
- We have spent the last 5 years developing an *in situ* forming, biodegradable polymer to provide two weeks of continuous retinal tamponade
- No positioning required, no induced refractive shift, no removal surgery
- Human pilot studies to begin Q3 2020

Current tamponade methods are a burden



Many promising alternatives have not advanced to the clinic


- Heavy silicone oils
- Magnetic oils
- Combination gas-oils
- Retinal sealants
- Intraocular capsules
- Hydrogels




Our Approach

- Find really smart polymer chemists
- Design a tamponade that meets design criteria
- Test it in animals for safety and efficacy
- Manufacture in compliance with Regulatory and Quality standards of the FDA
- Perform clinical trial

Our Approach



PYK-1105 (functionalized PVA-PEG hydrogel)



Physical Characteristics of PYK-1105

Parameter	Result
Time to injection (duration of minimal viscosity after mixing at 25°C)	10 min
Time to gel formation after injection (duration until viscous gel formation at 37 °C)	4 min
Time until degradation	11-14 days
Refractive index	1.3385
Transparency	>90% across visual spectrum

ISO 10993 Biocompatibility Data


- Evaluated in mice and rabbits using standard ISO protocols for evaluating biocompatibility of medical devices

Test	Result
Cytotoxicity	Non-cytotoxic
Intravitreal injection	Non-irritant
Acute systemic toxicity	Non-toxic
Material-mediated pyrogenicity	Non-pyrogenic
Endotoxin	< 0.2 EU/device
Implantation	No macroscopic reaction; minimal to no microscopic reaction
Genotoxicity	Non-mutagenic

Preclinical studies

- PYK-1105 formulation is well tolerated in a rabbit & mini-pig vitrectomy model
 - 25g vitrectomy/fluid-air exchange/PYK-1105 implantation
- Pre-GLP study in rabbit (4) & mini-pig (5)
 - One pig removed from study on POD3 due to post-op VH
- Manuscript accepted for publication in *JVRD*

Preclinical studies



Representative animals from
post-operative week 4

Preclinical studies

Rabbits

ID	Inflammation Score (conj redness/chemosis + AC reaction + iris)				Notes
	POD7	POD14	POD21	POD28	
PIK-35	3	3	0	0	POD7/14 score for conj injection
PIK-36	2.5	0	0	0	POD7 score for conj injection
PIK-37	1	0	0	0	POD7 score for conj injection
PIK-38	0	0	0	0	

	Vitreous Haze Score				
	POD7	POD15	POD22	POD28	
PIK-35	0	0	0	0	
PIK-36	0.5	0	0	0	
PIK-37	0.5	0.5	0	0	
PIK-38	0	0	0	0	

	Retina score (optic nerve + vasculature + retina + choroid)				
	POD7	POD15	POD22	POD28	
PIK-35	0	0	0	0	
PIK-36	0	0	0	0	
PIK-37	0	0	0	0	
PIK-38	0	0	0	0	

Pigs


ID	Inflammation Score (conj redness/chemosis + AC reaction + iris)				Notes
	POD7	POD14	POD21	POD28	
6-2	0	0	0	0	
9-2	0	0	0	0	
9-11	0	0	0	0	
4-2	0		0	0	no POD14 exam
89-1	0		0	0	no POD14 exam

	Vitreous Haze Score				
	POD7	POD14	POD21	POD28	
6-2	0	0.5	0	0	
9-2	0.5	0	0.5	0	
9-11	2	1	0.5	0.5	
4-2	0.5		0	0	no POD14 exam
89-1	1		1	1	no POD14 exam


	Retina score (optic nerve + vasculature + retina + choroid)				
	POD7	POD14	POD21	POD28	
6-2	0	0	0	0	
9-2	0	0	0	0	
9-11	1	0	0	0	
4-2	0		0	0	no POD14 exam
89-1	0		0	0	no POD14 exam

Clinical exam scores over first month
(Modified Draize/ McDonald-Shadduck scoring systems)

Preclinical studies




Representative animals from
Post-operative week 4



IOP remains normal
post-operatively



Normal histopath; normal ffERG





Rabbit H&E, post-operative week 4

Mini-pig H&E, post-operative week 4

Rabbit, pow4



Mini-Pig, pow4



Pilot human study

- 10 patient, first-in-human study to begin Q3 2020
- Multicenter, single cohort, open label study
- Enrollment restricted to patients with RRD who have limited visual potential
- Primary outcome is safety and tolerability
- Patients will be followed for 6 months post-operatively
- Interim study results expected Q1 2021

Team



James A. Stefater, MD, PhD
President & Cofounder
Vitreoretinal surgeon; Eye Health Services



Tomasz P. Stryjewski, MD
Chief Scientific Officer & Cofounder
Vitreoretinal surgeon; Tallman Eye
Co-Founder, Helio Vision Inc.



Sameer Sabir
Executive Chairman
Founder & CEO, Brixton
Founder & CEO, Seven Oaks
Founder & CEO, MoMelon

Leadership Team



Larry Roth
Product Development
20+ years experience in medical device development



Gordon Roberts
Quality
20+ years in medical device quality and regulatory



Maureen O'Connell
Regulatory
20+ years in ophthalmic device regulatory



Olivier Kagan
Project Management
Previously Director of Quality Systems at NSF



Lori Gilmartin, RN
Clinical Operations
20+ years in clinical study management

Advisors



Carl Awh, MD
Tennessee
Retina



Dean Elliott, MD
Mass Eye and Ear



Leo Kim, MD, PhD
Mass Eye and Ear



Darius Moshfeghi, MD
Stanford

Contact

Tony Stefater, MD, PhD
tony@pykustherapeutics.com

Tommy Stryjewski, MD
tommy@pykustherapeutics.com

Pykus Therapeutics Inc.
One Mifflin Place, Suite 320
Cambridge MA 02138